MI-5

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INFORMATION

PRODUCT

Mini incubator with 6 cavities for biological indicator - Box with 1 unit

PRESENTATION

Transparent lid.

CHARACTERISTICS

The MI-5 mini incubator consists of 6 cavities for simultaneous incubation of up to 6 biological indicators in the temperature range of 55 to 60°C. It has an exclusive system to break ampoules and a translucent protective lid providing more safety and hygiene during incubation. Its incubation chamber is made of aluminum. Automatic dual voltage, no need to switch the voltage key.

USAGE INDICATIONS

Indicated for the incubation of biological indicators between 55 and 60°C.

SPECIFICATIONS AND COMPOSITION

Incubation capacity: up to 6 indicators

Incubation temperature: 55 to 60°C

Incubation chamber composition: aluminum

Fairing composition: ABS Plastic

Voltage: Dual voltage

Power: 10 watts

Frequency: 50/60 Hz

STORAGE AND TRANSPORT

STORAGE AND WARNINGS

Store the product in a dry place, away from sunlight.

Never touch the incubation chamber. This area is heated and can cause burns.

TRANSPORT - RISK CLASSIFICATION

This product is not classified as dangerous according to the ANTT resolution n° 420 of February 12, 2004.

USER INSTRUCTIONS

Plug the power cord into the outlet (110-220V). At this time, the indicator LED will light up, indicating that the equipment is on.

After 15 minutes of heating, the MI-5 mini incubator is ready to be used. After this heating period, the self-contained biological indicators can be incubated.

Open the lid and insert the biological indicators. To carry out the microbiological monitoring, follow the manufacturer's instructions for the biological indicators.

APPLIED REGULATIONS

ANSI/AAMI/ISO 11138-1:2017, Sterilization of health care products Biological Indicators Part 1: General requirements.

EN 866-3:2000, Biological Systems for resting sterilizers and sterilization processes Part 3: particular systems for use moist heat sterilizers.